

Art-Nr.	Typ	Ges.Gew. in kg.	Eigengew. in kg.	L/B/H in mm ca.	Gebremst
040.020.367.007	Thermocase 1325/154 neu	1300	520	2520/1540/1500	Ja
040.020.367.008	Thermocase 1330/207 neu	1300	630	3090/2070/2030	Ja
040.025.367.013	Thermocase 1625/154 neu	1600	540	2520/1540/1500	Ja
040.025.367.019	Thermocase 1630/174 neu	1500	590	3090/1740/1930	Ja
040.025.367.016	Thermocase 1630/207 neu	1600	645	3090/2080/2030	Ja
040.025.367.018	Thermocase 1630/245 neu	1600	715	3090/2450/2230	Ja
040.025.367.020	Thermocase 1640/174 neu	1600	720	4090/1740/1930	Ja
040.025.367.021	Thermocase 1640/207 neu	1600	800	4090/2070/2030	Ja
040.035.367.026	Thermocase 2025/154 neu	2000	590	2520/1540/1500	Ja
040.035.367.027	Thermocase 2030/174 neu	2000	680	3090/1740/1930	Ja
040.035.367.031	Thermocase 2030/207 neu	2000	760	3090/2070/2030	Ja
040.035.367.034	Thermocase 2030/245 neu	2000	830	3090/2450/2230	Ja
040.035.367.028	Thermocase 2040/174 neu	2000	775	4090/1740/1930	Ja
040.035.367.032	Thermocase 2040/207 neu	2000	840	4090/2070/2030	Ja
040.035.367.035	Thermocase 2040/245 neu	2000	910	4090/2450/2230	Ja
040.045.367.023	Thermocase 2725/154 neu	2700	620	2520/1540/1500	Ja
040.045.367.064	Thermocase 2730/174 neu	2700	685	3090/1740/1930	Ja
040.045.367.068	Thermocase 2730/207 neu	2700	773	3090/2070/2030	Ja
040.045.367.073	Thermocase 2730/245 neu	2700	845	3090/2450/2230	Ja
040.045.367.065	Thermocase 2740/174 neu	2700	780	4090/1740/1930	Ja
040.045.367.072	Thermocase 2740/207 neu	2700	907	4090/2070/2030	Ja
040.045.367.062	Thermocase 2740/245 neu	2700	925	4090/2450/2230	Ja
040.050.367.030	Thermocase 3030/174 neu	3000	735	3090/1740/1930	Ja
040.050.367.035	Thermocase 3030/207 neu	3000	855	3090/2070/2030	Ja
040.050.367.040	Thermocase 3030/245 neu	3000	925	3090/2450/2230	Ja
040.050.367.031	Thermocase 3040/174 neu	3000	795	4090/1740/1930	Ja
040.050.367.036	Thermocase 3040/207 neu	3000	925	4090/2070/2030	Ja
040.050.367.041	Thermocase 3040/245 neu	3000	995	4090/2450/2230	Ja
040.055.367.030	Thermocase 3530/174 neu	3500	750	3090/1740/1930	Ja
040.055.367.034	Thermocase 3530/207 neu	3500	870	3090/2070/2030	Ja

040.055.367.039	Thermocase 3530/245 neu	3500	940	3090/2450/2230	Ja
040.055.367.031	Thermocase 3540/174 neu	3500	810	4090/1740/1930	Ja
040.055.367.035	Thermocase 3540/207 neu	3500	940	4090/2070/2030	Ja
040.055.367.040	Thermocase 3540/245 neu	3500	1010	4090/2450/2230	Ja

